

Accuracy and Stability of a Novel Subcutaneous Glucose Sensor

Jake Leach, John B. Welsh, Diana Le, Jonathan Hughes

Research & Development (JL, JH), Clinical Affairs (JBW), and Biometrics (DL), Dexcom, Inc., San Diego, CA, USA

Background and Aims

- CGM users value sensor performance characteristics such as accuracy, calibration frequency, and safety.
- We developed a new sensor and tested its accuracy with once-daily SMBG calibrations versus reference venous (YSI) values collected during clinic visits throughout 10 days of wear.

Methods

- Fifty subjects (49 with type 1 diabetes, mean±SD age 32.5±18.7 years) enrolled and 1 subject withdrew early.
- Subjects wore 1 sensor each, and used SMBG values for once-daily calibrations.
- CGM values were compared to YSI values collected on days 1, 2, 4, 7, and 10 of sensor wear.
- This analysis extends results in the abstract to include 994 paired YSI-CGM values from 49 subjects.

Results

- The overall percentage of sensor values within 15% of the reference value (or within 15 mg/dL [0.8 mmol/L] for reference values ≤100 mg/dL [≤5.6 mmol/L]) was **89.9% (95% CI, 87.9-91.7%)** and the analogous within-20% agreement rate was **96.1% (95% CI, 94.7-97.2%)** (Figure 1).
- The overall mean ARD was **8.1% (95% CI, 7.6-8.5%)** (Figure 2). The median ARD was **6.8%**.
- Accuracy by all 3 metrics was best on Day 2.
- Surveillance error grid analysis (Figure 3) showed that almost all of the YSI-CGM paired points were in the "No Risk" or "Slight Risk" zones.
- The outlier point (YSI = 183 mg/dL, CGM = 70 mg/dL) would not prompt an insulin bolus if the CGM datum was used for an insulin dosing decision.




Figure 1. Within 15% and within-20% agreement rates




Figure 2. Mean absolute relative differences




Figure 3. Surveillance error grid and risk categories

Summary

- The new system is accurate throughout its 10-day functional life.
- Once-daily calibrations may improve the user experience by reducing the need for SMBG.
- The surveillance error grid analysis supports the new system's potential utility for supporting diabetes management decisions and in closed-loop insulin delivery systems.